

## Acetabular Reamers



The acetabular reamer is used in total hip arthroplasty or replacement to prepare the acetabulum to receive the replacement liner. Greatbatch Medical's cross design was the first to allow for a universal handle for all reamer sizes.

### INNOVATIVE DESIGN FOR EFFICIENT PROCEDURES

- Patented cross design allows for various handle and reamer combinations.
- Enables precision spherical reaming.
- Easy collection of reamed bone stock.

### A CUT YOU CAN COUNT ON

- Secure and easy to maneuver, providing a more accurate cut.
- Available in standard tooth, contour tooth, full hemisphere and MIS versions.
- Hemispherical reamers provided, cut in a range of 36 to 80 mm.

---

**Greatbatch Medical Headquarters**  
10000 Wehrle Drive  
Clarence, NY 14031  
tel 716-759-5600  
fax 716-759-2562

**US Customer Service**  
4532 Park 30 Drive  
Columbia, City, IN 46725  
tel 260-244-6300  
fax 260-244-6345

**Europe & Asia-Pacific Customer Service**  
L'Echelette 7  
2534 Orvin  
Switzerland  
tel +41 32 358 01 11  
fax +41 32 358 01 00

[orthopaedicsna@greatbatchmedical.com](mailto:orthopaedicsna@greatbatchmedical.com)

[orthopaedicseurope@greatbatchmedical.com](mailto:orthopaedicseurope@greatbatchmedical.com)

[www.greatbatchmedical.com](http://www.greatbatchmedical.com)